

## Congratulations on choosing the Home & Yard System as your personal emergency response system!

The Home & Yard System is easy to install. Simply follow these step-by-step instructions to get started using your new personal emergency response system.

1. Unpack the Home & Yard System from its box and check that you have all the required parts before you begin installation of the system. We recommend that you keep the system's box in a safe place in case you ever need to ship it back for repair or return.

*\*\*\*Please fill out all of the enclosed paperwork and set it aside to be mailed in to the appropriate address provided.*

2. Select a room in your home where you would like to set up your personal emergency response system. We encourage you to avoid placing the system next to televisions, radios, air conditioners, and other noisy devices to ensure you will be heard easily by the Galaxy Medical Alert Response Center in an emergency. The room must have a telephone jack and an electrical outlet that can be used for the system (it is important that the system is connected to an electrical outlet that is not controlled by a wall/light switch). Connect the provided AC Adaptor to the DC 7.5 Volt jack on the back of the system console. Plug the AC Adaptor into an electrical outlet that is **not** controlled by a wall/light switch.
3. Turn the ON/OFF switch on the underside of the system console to the ON position.
4. Unplug the telephone line from the back of your telephone and plug it into the Line In jack (A) on the system console; the phone cord coming from the wall should now be connected to the console. Plug one end of the telephone cord received with your console into the To House Phone jack (B) on your console and the other into the back of your telephone.
5. It is now time to test your Home & Yard System. To test the system, press the blue button on your personal emergency response **pendant** and hold for 4 seconds. After a series of beeps and a dialing sound, you may experience approximately 30-60 seconds of silence followed by a response from the Galaxy Medical Alert Response Center asking if you are in need of assistance. Simply inform them that you are testing your new personal emergency response system. The Galaxy Medical Alert Response Center may proceed to ask you some questions to further set up your account. Set up of the Home & Yard System is complete once the Galaxy Medical Alert Response Center greets you as a new client and disconnects from the call.

**Note:** *If the Galaxy Medical Alert Response Center does not respond within 2 minutes of the first test call, please call 1-855-905-7414 for technical assistance.*

6. Please return the completed paperwork in the stamped envelope that was provided with the Home & Yard System. If you require any assistance with filling out the form(s), please call us at 1-855-905-7414.

We recommend that you test your Home & Yard System monthly. To test your system, simply repeat step 5 above once a month. You do NOT need to call us before testing. If at anytime you are unable to make a successful test call to the Galaxy Medical Alert Response Center, please call 1-855-905-7414 for assistance.

