

Quick Reference

Turn ON: Press the Call Button and hold until the Battery Indicator flashes.

Turn OFF: Press the Call Button twice quickly; hold the second press for 7 seconds or until the system plays a message confirming your request to turn off. Press the Call Button again to confirm.

Place an emergency call:

1. Press and hold the Call Button.
2. Listen for the alert tone, followed by the audio message: "Calling your response team now. Please hold the device and position the microphone near your mouth to complete your call."
3. Your response team will answer and notify emergency services for you.

Battery & Cellular Indicators

Battery Indicator around Call Button shows the level of battery charge.

Battery Indicator	Status*	Action Required
Green	Good	Ready for use
Amber	Moderate	No action needed
Red	Low	Charge soon

*Battery is designed to hold a charge for 36 hours but may be reduced as a result of user activity, cellular coverage, talk time, device-specific settings, and device life.

Cellular Indicator on top of the system shows cellular signal strength.

Cellular Indicator	Signal Strength*
Green	Strong
Amber	Moderate
Red	No Signal

*A cellular signal may not be available in all areas.

About The Home & Away System with GPS



1. **Cellular Indicator** shows cellular signal strength.
2. **Battery Indicator** shows battery status.
3. **Call Button** initiates a two-way call to response team.



Phone: 1-855-905-7414
info@galaxymedicalalert.com

www.galaxymedicalalert.com

OUR LOCATIONS

30 Rayette Road
Concord, Ontario
Canada L4K 2G3

361 Joffre Place
Windsor, Ontario
Canada N8X 2T1

Home & Away System with GPS

Quick Start Guide



Enjoy New Freedom

This Quick Start Guide will help you set up and activate the Home & Away System with GPS. **Please read this entire guide before using the Home & Away System with GPS to ensure it is set up correctly and ready for use.**

What's in the Box



Step 1: Activate the System

1. Insert the Power Cord into the Charger. Plug the large end of the Power Cord into a standard wall outlet.
2. Place the system in the Charger.
3. The system will automatically turn on and begin activation.
4. Follow the voice prompts until directed to connect to your response team.



Step 2: Call to Set Up

1. Press and hold the Call Button.
2. Your response team will answer shortly to customize your system and complete set up.
3. After the call, charge the system in the Charger for 4 hours. Then it is ready for use.

IMPORTANT:

The unit is **not ready for use** until you:

- Charge system in the charger for 4 hours
- Complete a call with the response team



Step 3: Wearing the System

The Home & Away System with GPS can be worn on a belt using the Belt Clip, or as a pendant around your neck using the Lanyard/Clip. If worn with the Lanyard, the system should be worn under clothing to prevent excess movement of the device. To attach the Belt Clip or Lanyard/Clip:

1. Align the notched side of Clip with the notched area on the back.
2. Press inward until the Clip notches are nested within the system's notches.
3. Firmly push the Clip toward the top of the system. It is secure when you hear a click or you can no longer see the notches on the Clip.



Step 4: Review User Guide

The User Guide provides key information about the Home & Away System with GPS, including Important Safety Information. Please review it carefully and keep it nearby for reference.

If you have questions, call your response team.

By activating your Home & Away System with GPS you agree to the Terms of Use and understand the Important Safety Information in the the Home & Away System with GPS User Guide.

